









# POSITIVE AIRWAY PRESSURE SYSTEMS

Company		3B Medical Inc	Drive DeVilbiss Healthcare	Fisher & Paykel Healthcare	Human Design Medical, LLC
<b>Product</b>		 Luna CPAP & Auto-CPAP	 IntelliPAP AutoAdjust	 ICON+ Series	 The Z1 Auto
<b>Website</b>		www.3Bproducts.com	www.drivemedical.com	www.fphcare.com	www.smallcpap.com
<b>Warranty (years)</b>		2	3	2	3
<b>Operating Modes</b>	<b>CPAP</b>	Yes	Yes	Yes	Yes
	<b>APAP</b>	Yes	Yes	Yes	Yes
	<b>Bilevel</b>	No	No	No	No
	<b>Auto-adjusting Bilevel</b>	No	No	No	No
<b>Pressure Relief</b>		Yes	Yes	Yes	Yes
<b>Dimensions (cm)</b>		17 x 19 x 11	10.7 x 16.5 x 17.5	16 x 17 x 22	16.5 x 8.3 x 5.1
<b>Weight (kg)</b>		1.5	1.22	2.2	0.28
<b>Noise Level (dBA)</b>		28	26	<30	26
<b>Operating Pressure Range (cm H<sub>2</sub>O)</b>		0-20	3-25	4-20	4-20
<b>Pressure Adjustment Increment (cm H<sub>2</sub>O)</b>		0.5	0.5	0.5	0.5
<b>Adjustable Starting Ramp Pressure (cm H<sub>2</sub>O)</b>		4	3-25	4-20 (SensAwake)	4.0-19.5
<b>Ramp Time (minutes)</b>		0-60	0-45	Automatic (SensAwake) or 20 minutes	0-45
<b>Data Recorded</b>		Usage, leak, AHI, pressure, 95th percentile pressure, snore, central apnea index, compliance percentage	Usage, AHI, leak, pressure, exhale puffing, non-responding events, 90th and 95th percentile pressures, snore, hypopnea, adherence score, oximetry (optional)	Hours used, compliance rate, AHI, leak, AI, HI, flow limitations, detailed night-by-night data, 90th percentile pressure, SensAwake detections, mask on/off	Usage, leak, AH count, AHI (index)
<b>Features</b>	<b>Auto Start/Stop</b>	Yes	Yes	No	Yes
	<b>Leak Compensation</b>	Yes	Yes	Yes	Yes
	<b>Mask On/Off Alert</b>	Yes	No	Yes	Yes
	<b>Altitude Compensation</b>	Yes	Yes	Yes	Yes
	<b>Heated Humidification</b>	Yes	Yes	Yes	No
	<b>Remote Programming</b>	No	Yes	Yes	Yes
<b>Patient Access to Own Data Without Physician Involvement</b>		Yes, via free smartphone app, iCodeConnect Quick Report, and on screen data.	Yes, via 3 options. 1. Entering the SmartCode from the device display at www.intelliPAP.com. 2. Setting up an account at www.IntelliPAPSmartCode.com. 3. Using the QuickView menu on the device directly.	Yes, via on-screen patient data menu.	Yes, via HDM's Nitelog App, available for free on iOS and Android.
<b>Additional Features</b>		WiFi wireless data and QR code scanning for smartphone app.	Made in the USA, built in Somerset, Pa, and has less than a 1% failure rate over 5 years; optional 2-year extended warranty to cover the patient for a full 5 years. The IntelliPAP platform offers a full line of devices: Standard, Standard Plus, AutoAdjust, Bilevel S, and AutoBilevel.	Aesthetically pleasing; simple mode for patients and clinicians; ThermoSmart humidification; SensAwake responsive pressure relief; InfoSmart Web cloud-based software; flexible data methods, eg, InfoUSB, InfoGSM modems, and InfoView; reliability. The F&P ICON+ family includes Novo, Premo, and Auto.	Integrated overnight battery.

Company		Human Design Medical, LLC	Philips Respironics	ResMed	Somnetics International Inc
<b>Product</b>		 The Z1 CPAP	 DreamStation	 AirSense 10 AutoSet	 Transcend Auto
<b>Website</b>		www.smallcpap.com	www.philips.com/dreamstation	www.resmed.com/us/en/consumer/products/devices/airsense-10-auto-set.html	www.mytranscend.com
<b>Warranty (years)</b>		3	2	2	3
<b>Operating Modes</b>	CPAP	Yes	Yes	Yes	Yes
	APAP	No	Yes	Yes	Yes
	Bilevel	No	Yes	No	No
	Auto-adjusting Bilevel	No	Yes	No	No
<b>Pressure Relief</b>		Yes	Yes	Yes	Yes
<b>Dimensions (cm)</b>		16.5 x 8.3 x 5.1	15.7 x 19.3 x 8.4	11.6 x 25.5 x 15.0	17.6 x 10.2 x 8.7
<b>Weight (kg)</b>		0.28	1.33 (without humidifier); 1.98 kg (with humidifier)	1.25	0.43
<b>Noise Level (dBA)</b>		26	25.8 (device); 27.9 (device with humidifier)	25	26.6
<b>Operating Pressure Range (cm H<sub>2</sub>O)</b>		4-20	4-20	4-20	4-20
<b>Pressure Adjustment Increment (cm H<sub>2</sub>O)</b>		0.5	0.5	0.2	0.1
<b>Adjustable Starting Ramp Pressure (cm H<sub>2</sub>O)</b>		4.0-19.5	4	4-20	4-20
<b>Ramp Time (minutes)</b>		0-45	0-45 (5 minute increments)	0-45 min or auto (begins at sleep onset)	0-45, adjustable
<b>Data Recorded</b>		Usage, leak, AH count, AHI (index)	Encore data capture and reporting, use patterns, advanced event detection,* flow waveform reports. *This can detect, respond to (some models), and log: AHI, RERA, snore, flow limitation, leak, vibratory snore, obstructive airway apnea, large leak, clear airway apnea, periodic breathing, hypopnea.	Hours used, mode of therapy, apnea (central, obstructive) and hypopnea (AHI), leak, CSR, RERA, pressure, EPR level	Average hours/night, median hours/night, therapy session (4, 6, and 8 hours), therapy settings, time-weighted average pressure, 90th and 95th percentile pressure, AHI index, apnea index, hypopnea index, leak, nightly detailed compliance including: pressure trend, leak rate, apneas, hypopneas
<b>Features</b>	Auto Start/Stop	No	Yes	Yes	No
	Leak Compensation	Yes	Yes	Yes	Yes
	Mask On/Off Alert	Yes	No	No	No
	Altitude Compensation	Yes	Yes	Yes	Yes
	Heated Humidification	No	Yes	Yes	Yes
	Remote Programming	Yes	Yes	Yes	No
<b>Patient Access to Own Data Without Physician Involvement</b>		Yes, via HDM Data Viewer Software.	Yes, via the DreamMapper patient self-management tool.	Yes, wirelessly via myAir and on the device screen "Sleep Report."	Yes, download the free software, then connect device via provided mini USB cord. Or use the TranSync cloud system, which also requires a mini USB connection.
<b>Additional Features</b>		Integrated overnight battery.	Opti-Start feature, EZ-Start Auto-Trial, C-Flex, C-Flex+, A-Flex, Humidity Control; Heated Tube; Advanced Event Detection; Performance Check EncoreAnywhere; DreamMapper; cellular or WiFi modem accessories, SD card; optional oximetry accessory.	Integrated humidifier, built-in wireless data communications, Advanced Event Detection, AutoRamp, Climate Control, SmartStart, optional integrated oximetry, Remote Assist, AirView, myAir. AirSense 10 is also available in CPAP, Elite, and AutoSet for Her models.	Multiple portable battery options—rechargeable with AC, DC, and solar power—to make it easy for patients to use CPAP therapy anywhere.