

Company	3B Medical Inc	AG Industries LLC	Airway Management Inc	CareFusion	Circadiance
<b>Product</b>	 iVolve MiniNasal	 Nonny	 myTAP PAP Nasal Pillow Mask	 Advantage Series Full Face Mask	 SleepWeaver ADVANCE Soft Cloth Nasal CPAP Mask
<b>Website</b>	www.3bproducts.com	www.agindustries.com	www.mytappap.com	www.carefusion.com	www.circadiance.com
<b>Warranty (months)</b>	3	1	3 (patient compliance guarantee)	1	3
<b>Available Sizes</b>	S, M, L	S, M, L, XL (small adult)	S, M, L nasal pillows included	S, M, L, XL	S, one size fits most
<b>Latex-free</b>	X	X	X	X	X
<b>Intentional Leak @ 10 cm H<sub>2</sub>O (l/min)</b>	29	30	30	Info not available	25
<b>Fitting Support Solutions</b>	Comes with instructional booklet with fitting instructions and troubleshooting guide.	Fitting instructions included in mask kit.	Patient and clinician fitting guide included; fitting video on myTAP-PAP.com and DVD available.	Instructions for use included with purchase.	www.circadiance.com; DVD is included with the mask.
<b>Pressure Range (cm H<sub>2</sub>O)</b>	4-30 (recommended)	4 – 20	4-18	3-20	4-20
<b>Cleaning Instructions</b>	Hand wash by gently scrubbing in warm water using mild soap. Use a soft-bristled brush to clean the port cover. Rinse all components well with water and allow to air dry out of direct sunlight.	Wipe off cushion daily using wet cloth or Mask Wipe. Once a week, disassemble components and immerse in warm soapy water. Rinse well and dry before use.	After each use, rinse the mouth-piece in cold water and wipe the pillows' outer portion with a damp cloth. Once a week, disassemble the mask and soak parts for no longer than 10 minutes in cold water and mild liquid soap. Remove and rinse in cold water then allow the parts to air dry.	Clean with mild detergent. It does not require disassembly.	Clean by hand washing with mild nonabrasive soap in warm water. Rinse well. Hang to dry.
<b>Headgear Description</b>	Quick-release headgear with clips that simplify attachment and release.	5-point adjustable straps to ensure custom fit. Easy to release and apply quick clips.	Headgear is available for the myTAP PAP Nasal Pillow Mask. However, the mask does stay securely in place without straps.	Adjustable headgear.	Headgear is soft, comfortable, and features a 3-point adjustable strap system.
<b>Interface Description</b>	Composed of only 2 parts, the mask is extremely easy to fit. Its lightweight and small profile makes it a perfect mask for patients looking for a lighter and less obtrusive mask.	Mask Kit comes with 2 size cushions, which can be interchanged onto mask frame to accommodate developing faces. Soft medical grade silicone, replaceable mask cushions provide a comfortable fit. Easily adjustable mask frame allows 3 mask positions for custom fit to facial features.	The myTAP PAP Nasal Pillow Mask is a strapless nasal pillow mask that is intended for individuals who have been diagnosed as requiring CPAP or bi-level ventilator treatment. The myTAP PAP is intended for single patient adult use in the home, hospital, or other clinical setting.	Seal around the nose and mouth with triple layer seal design.	The SleepWeaver uses the properties of a balloon. When a balloon is inflated, all points inside the balloon are at the same pressure. This allows the surface areas of the mask to conform to the face.
<b>Additional Features</b>	The iVolve MiniNasal has quiet venting and very soft silicone cushioning to improve patient comfort.	Soft flexible tubing attached to mask to allow less torque for restless sleepers when changing positions while sleeping.	The myTAP PAP Nasal Pillow Mask encourages patient compliance by solving the problems commonly associated with traditional CPAP Masks—from awkward straps to disruptive air leaks. The unique and innovative design of the myTAP PAP NPM delivers unparalleled performance and is backed by a Compliance Guarantee Program.	Elbow swivel 360 degrees with anti-aphyxia valve, built-in vent openings.	SleepWeaver lets you fit your patients quickly and comfortably to attain a reliable reading during titration.

Information for this guide based on data submitted by product manufacturers. Sleep Review strives for accuracy in all data but cannot be held responsible for claims made by manufacturers.

Fisher & Paykel Healthcare	Hans Rudolph Inc		InnoMed Technologies	Philips Respironics		ResMed
						
Simplus Full Face Mask	V2 Full Face CPAP Mask	QUEST HOME CARE Full Face Mask	Aloha Nasal Pillow	DreamWear Under the Nose Nasal Mask	Amara View Minimal Contact Full Face Mask	AirFit F10
<a href="http://www.fphcare.com/simplus">www.fphcare.com/simplus</a>	<a href="http://www.rudolphkc.com">www.rudolphkc.com</a>	<a href="http://www.rudolphkc.com">www.rudolphkc.com</a>	<a href="http://www.innomedinc.com">www.innomedinc.com</a>	<a href="http://philips.com/dreamwear">philips.com/dreamwear</a>	<a href="http://philips.com/amaraview">philips.com/amaraview</a>	<a href="http://www.resmed.com/us/en/consumer/products/masks/airfit-f10.html">www.resmed.com/us/en/consumer/products/masks/airfit-f10.html</a>
3	6 (includes size guarantee)	6	3	3	3	3
S, M, L	P, XS, S, M, L	P, XS, S, M, L	S, M, L	S, M, L, MW	S, M, L	XS, S, M, L
X	X	X	X	Not made with natural rubber latex	Not made with natural rubber latex	X
35	N/A	No, pressure tested	28	26.6	26.2	37
Fitting gauge included; videos at <a href="http://www.youtube.com/fphcare">www.youtube.com/fphcare</a> .	Call or e-mail HRI for pictured fitting instructions.	Mask sizing gauge, fitting instructions included; videos at <a href="http://www.rudolphkc.com">www.rudolphkc.com</a> .	Comes with all sizes. Includes step-by-step fitting instructions. For more information, call (800) 200-9842.	Quick start guide attached to each mask, full instructions in OIS and support tools in DreamMapper.	Full instructions in OIS and support tools in DreamMapper.	See <a href="http://www.resmed.com/us/en/consumer/support/masks/how-to-fit-your-mask.html">www.resmed.com/us/en/consumer/support/masks/how-to-fit-your-mask.html</a> .
4-25	3-40	3-40	4-20	4-20	4-30	3-40
Minimal parts to clean; soap in lukewarm water; high level disinfection available.	Mask, swivel dishwasher safe. Can hand-clean headgear. V2: steam autoclaved, liquid chemical, hot water pasteurized in a hospital setting.	Soak components in warm water with mild detergent, sponge or soft brush, rinse in warm water, air dry, or use lint-free cloth.	Wash daily with warm and mild soapy water. Rinse well. Air dry out of direct sunlight.	Hand wash in warm water with a mild liquid dish washing detergent. Rinse thoroughly. Air dry completely before use. Lay the fabric parts flat or line dry.	Hand wash in warm water with a mild liquid dish washing detergent. Rinse thoroughly. Air dry completely before use. Make sure the entire mask is dry before use. Lay the fabric parts flat or line dry.	Gently hand wash with mild soap; avoid harsh cleaning products; allow to dry away from direct sunlight.
ErgoForm breathable headgear incorporates stretch and non-stretch panels. It holds its shape and allows stretch where required for optimal usability and performance. Embedded Velcro tabs with easy-hook clip hooks enable fast fitting and removal. Crown strap for stability.	Simple strap foam headgear with side adjustments and top of the head tri-glide adjustment.	Light gray-on-gray color; soft foam and nylon; mounts with simple 4 clips snapped into place; secure comfortable fit.	Comes standard with Comfort Wraps. Features easy push on/pull out reinforced headgear for connecting to Pillow Reservoir.	All fabric non-slip headgear.	All fabric headgear.	The seamless SoftEdge headgear is gentle on the skin and features premium rolled-edge fabric for reduced facial markings and easy fit. The SoftEdge headgear requires minimal tension to achieve a comfortable fit.
The F&P Simplus incorporates three key components: the RollFit Seal, ErgoForm Headgear, and Easy Frame, all designed to work in harmony. In combination, these components offer the comfort, seal, and easy use to simply fit and perform.	One piece of molded silicone. No painful forehead support. No hard plastic frame. Goes below the chin to the top of the nose.	Fits into the face securely for leak-free fit, won't ride up face due to chin cup, no hard frame so no skin breakdown, 5 sizes, glasses allowed while wearing.	The Aloha is an extremely lightweight and quiet nasal pillow. Offers an incredibly stable and comfortable seal. Has Floating Ball and Socket Elbow. Pillow angle can be optionally secured over-the-head.	DreamWear's unique design offers many of the benefits of nasal and pillows masks to allow patients to have the best of both mask types. Allowing more freedom of movement and providing exceptional comfort, DreamWear lets patients feel like they are not wearing a mask at all.	Amara View's incredible under-the-nose design makes it the only full face mask that won't cause red marks on the bridge of your nose when compared to the market leaders. Amara View covers less of a patient's face than all leading full face masks.*	The AirFit F10's four-part design makes it easy to take apart, clean, and reassemble. It is quiet, lightweight, and provides a clear line of sight. The dual wall Spring Air cushion enhances comfort at the nasal bridge and maintains a high performance seal.
Integrated diffuser that minimizes noise and draft for the bed partner, lightweight frame with easy clip, ball and socket swivel connector.	Detachable swivel port, airtight seal, new headgear.	Low cost home care mask allows all patients to afford to wear the high quality Hans CPAP/BiLevel mask. Distributor pricing available.	Innovative Arced-Track technology allows the angle of the Pillows Reservoir to be adjusted to a more customized, natural fit.	Tubing attachment located at the top of the head to offer freedom of movement during the night. Prevents red marks, discomfort, or irritation on the nostrils or on the nose bridge (disclaimer: mask does not directly contact the bridge of the nose or nostrils). Soft, flexible silicone frame and fabric wraps provide a soft touch on the patient's cheek.	Innovative design prevents red marks, discomfort, or irritation on the nose bridge. The smallest and lightest† of all leading full face masks. Widest field of vision of all leading full face masks.* Enables wearing glasses, reading, and watching TV. Easy to use with one modular frame for all cushion sizes. No bulky cushion and frame in front of your face and eyes. *Internal testing; data on file † Without quick release tube; data on file	AirFit F10 full face mask is comfortable, reliable, and easy to use. It has a sleek design, and is the lightest full face mask on the market. A high performance seal delivers effective therapy.